GLIADIN AND CELIAC DISEASE

RESULT/ERGEBNIS/RÉSULTAT/RISULTATO	PRIMA[®]
<input type="checkbox"/> P	HOME TEST
<input type="checkbox"/> N	
Cut-off: 20 U/ml	
<i>date / Datum / date / data</i>	
CE 0483	

Self-test to check for the presence of anti deamidated gliadin (anti-DGP) antibodies, IgA and IgG, in the blood, for celiac disease screening.

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LITERATURANGAßEN

1. PD. Mooney, S.H.Wong, A.Johnston, M.Kuren, A.Averinos and D.Sanders "Increased detection of celiac disease with measurement of Deamidated Gliadin Peptide Antibody before endoscopy" *Cl.*,Gastroenterology and Hepatology 2015;13:1278-1284

2. N. R. Lewis & B. B. Scott "Meta-analysis: deamidated gliadin peptide antibody and tissue transglutaminase antibody compared as screening tests for celiac disease" *Alimentary Pharmacology & Therapeutics*

3. A. Lamm, P. Arkoski, S. Simell, T. Kinnunen, V. Simell, S Paavanen-Huhtala, A. Hinkkanen, R. Veijola, M. Knip, J. Toppin, O. Vaarala, O. Simell, and J. Ilonen "Antibodies to Deamidated Gliadin Peptide in Diagnosis of Celiac Disease in Children" *Gastroenterology* - JPCN 2015;60: 626–631

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PRECAUTIONS

1. Read these instructions for use carefully before performing the test. The Test is reliable only if all the instructions are followed correctly.

2. Keep the Test out of the reach of children.

3. Do not use the Test after the expiry date or if the package has been damaged.

4. Follow the procedure exactly, using only the specified quantities of blood and diluent.

5. Store the Test components at a temperature between +4°C and +30°C. Do not freeze.

6. Use the test and lancing device once only.

7. The test is for external use only. DO NOT SWALLOW.

8. In vitro diagnostic device for people who take anti-coagulant medications (blood thinners) or people suffering from haemophilia.

10.After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.

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